

# BROMINE FLOATER

## Instructions For Use

*A dispenser for those who use bromine tablets as their primary sanitizer.*

### **IMPORTANT! Follow All Directions and Cautions on Chemical Bottle.**

1. Remove the cap on the large end of the floater by turning the cap counter-clockwise one quarter of a turn and lifting off.
2. Using caution against spillage, fill floater with bromine tablets. Replace the lid and turn clockwise one quarter turn until locked.
3. Open dispenser adjustment on bottom of floater to expose

approximately half the holes. Set the floater in the spa. After 24 hours, test the water to determine the bromine level. Test bromine levels daily and adjust the dispenser opening as needed to maintain the recommended level of bromine.

Remove the floater when using spa.